



SOLIS
BIODYNE

Product Catalogue 2019

Ambient Temperature Stable
PCR and qPCR Reagents





Greetings from the Solis BioDyne Team

Dear Partners,

The past year has been full of exciting developments for Solis BioDyne. Quality has always been the core value of our work. To ensure we match the high quality management requirements of our diagnostic sector partners, we implemented the ISO 13485 standard.

We have continued developing our reverse transcription portfolio and have released the first products in our one-step RT-PCR range, which you will already find in this catalogue.

Inspired by questions and discussions with our customers we will keep improving and hope to be the partner of choice in your molecular biology applications.

Best regards,

Kadri Artma
Chief Executive Officer

About us

Solis BioDyne manufactures high quality molecular biology reagents since 1995. The company is based in Tartu, Estonia, an international academic city with a growing biotechnology sector.

Our young and professional team is dedicated to provide:

- **high quality products**
- **cost effective solutions**
- **quick and personal service**

Our product line includes DNA polymerases, master mixes for PCR/qPCR, reverse transcriptases and other reagents – all stable at ambient temperature. Solis BioDyne Quality Management System is certified to ISO 13485:2016 and ISO 9001:2015 standards.



Table of Contents

Stability TAG Technology	2
qPCR Mix Specifications & Mix Compatibility Table	4
Dye-based qPCR Mixes	7
Probe-based qPCR Mixes	11
Product selection Guide: Standard PCR Enzymes and Master Mixes	14
Hot-start PCR Enzyme and Mixes	15
Standard PCR Enzyme and Mixes	20
NEW! Product selection Guide: One-step RT-PCR Kits	22
Probe-based one-step RT-qPCR Kits	23
Product selection guide: cDNA synthesis	25
Reverse transcriptases and cDNA synthesis kits	26
Other enzymes and Reagents: dNTP Mix and Set, DNA Ladders, DNA Loading Dye Buffers etc.	29
Ordering	33
Distributors	35
Product List	36



Teaduspargi 9, 50411 Tartu, Estonia
T: +372 7409 960
F: +372 7402 079
E: info@solisbiodyne.com
Skype: solis.biodyne
W: solisbiodyne.com


VAT No: EE100587614
Reg. No: 10242922

Bank details: Swedbank AS

IBAN code: EE692200221005142234, SWIFT/BIC.: HABAE2X
 Bank address: Liivalaia 8, 15040 Tallinn, Estonia

Polypeptide Stabilization Technology: Stability TAG


3D MODEL OF FIREPOL® DNA POLYMERASE INCLUDING STABILITY TAG



■ Taq DNA Polymerase ■ Stability TAG

Our technology

We have genetically modified our enzymes to improve their long term stability and enhance stability at ambient temperatures. The addition of the Stability TAG (EU Patent EP2501716 and US Patent No 9,321,999) ensures that our enzymes are fully active even after keeping them at room temperature for a month!



■ 5x HOT FIREPol® Probe Universal qPCR Mix -20°C ■ 5x HOT FIREPol® Probe Universal qPCR Mix room temperature

Our qPCR mix 5x HOT FIREPol® Probe Universal qPCR Mix shows no loss of activity after incubation at room temperature for 1 month!


Ice-free reaction set-up

Our temperature stable enzymes allow you to change the way you work with enzymes that are usually extremely sensitive to temperature changes.


- **more space** on your work station and in your freezer
- **convenient** working conditions
- **less energy used** for making ice or keeping bench-top coolers cold



REACTION SET-UP ON ICE



REACTION SET-UP AT ROOM TEMPERATURE



Worry-free storage

Our reagents remain fully active even after a power outage has damaged everything else in your freezer, after someone has forgotten the reagents on the table overnight, or if there have been delays in the customs during shipment.



Solis BioDyne products are excellent products and meet the following criteria:

- **quality:** they give good results and a sharp signal without any background noise.
- **cost:** they are inexpensive, the order process is fast and the fact they are storable/shipped at room temperature reduces the cost of transportation.
- **convenience:** no need to order separately the PCR products as everything is inside the Mix which reduces pipetting errors. It is not necessary to work in ice which helps to save energy and avoid the deterioration of product due to power cuts or forgetting the product at room temperature.




ABDUL AZIZ WANE

Biology and Sanitary Engineer

Institut Pasteur Dakar, Senegal

No need to pay
for shipping the ice




Ice-free shipping

Our Stability TAG technology means we can ship your order without dry ice and large insulation boxes. Which is better for the environment and your budget:

- **less material** used for packaging
- **lower package weight**
- **less fuel** used for transportation
- **lower shipping charges** for you

Find out how much you can save when you
don't have to pay for shipping the dry ice



Ca 4x savings

Ca 3x savings

Ca 2x savings

qPCR Mixes: Specifications

5x concentration

All Solis BioDyne PCR and qPCR mixes are produced in 5x concentrations to improve stability. The higher buffer concentration also leaves 2.5 times more room for your template and primers compared with 2x mixes. This is highly advantageous when working with low-concentration DNA samples or performing multiplex assays.

EvaGreen® and SolisGreen™

Solis BioDyne dye-based qPCR mixes contain next generation DNA binding dyes EvaGreen® or SolisGreen™. The fluorescence spectra of these dyes are similar to the more widely used SYBR® Green I and are compatible with all major real-time cyclers. Our next generation dyes have:

- high **fluorescence** level
- high **sensitivity** for detecting low template concentrations
- high **stability** for room temperature storage

Light-protective packaging

Fluorescent dyes EvaGreen®, SolisGreen™ and ROX are sensitive to degradation by light. Solis BioDyne qPCR mixes that contain these dyes are supplied in special dark vials to minimize light exposure during transportation and reaction set-up.



COMPARISON OF OUR 5X QPCR MIX TO A STANDARD 2X QPCR MIX

	Vial size	Reactions (20 µl final volume)
2x Standard qPCR Mix	1 ml	100 rxn
5x Solis BioDyne qPCR Mix	1 ml	250 rxn

» Get 2.5x more reactions done with same volume of qPCR Mix


Quality Control

To ensure the highest quality and lot-to-lot consistency, strict quality control procedures are followed and documented in all stages of production. Certificates of Analysis are available for each lot upon request. Some tests that are performed with each lot of qPCR mixes:

- functional qPCR test with different templates and primer pairs
- amplification efficiency $\geq 10^5$ fold
- no endodeoxyribonuclease, exodeoxyribonuclease or self-priming activities
- comparison with previous lots

HIGHLY COMPETITIVE

Four tenfold dilutions of the PPIA gene from human cDNA were amplified using 5x HOT FIREPol® SolisGreen™ qPCR Mix (red) and qPCR mixes from Competitor A (green) and Competitor B (blue). The amount of DNA per reaction ranges from 0.1 ng to 0.1 pg. The reactions were performed on an Applied BioSystems™ QuantStudio™ 6.



qPCR Mix Compatibility Table: Dye-based qPCR Mixes

qPCR Platforms	5x HOT FIREPol® SolisGreen™ qPCR Mix	5x HOT FIREPol® EvaGreen® qPCR Supremix	5x HOT FIREPol® EvaGreen® qPCR Mix Plus (ROX)	5x HOT FIREPol® EvaGreen® qPCR Mix Plus (no ROX)	5x HOT FIREPol® EvaGreen® qPCR Mix Plus (Capillary)	5x HOT FIREPol® EvaGreen® HRM Mix (ROX)	5x HOT FIREPol® EvaGreen® HRM Mix (no ROX)
Applied Biosystems: 5700, 7000, 7300, 7700, 7900HT, StepOne™, StepOnePlus™		●	●			●	
Applied Biosystems: 7500, ViiA™7, QuantStudio™ 3, 5, 6 Flex, 7 Flex, 12K Flex	●	●	●			●	
Agilent/Stratagene: Mx3000P™, Mx3005P™, Mx4000™	●	●	●				
BioRad: CFX96™, CFX384™	●	●		●			●
BioRad: iQ™5, MyiQ™, Chromo4™, Opticon®2; MiniOpticon®		●		●			●
Bio Molecular Systems: Mic	●	●		●			●
Eppendorf: Mastercycler® ep Realplex	●	●		●			
Qiagen: Rotor-gene® 3000, Rotor-gene® 6000, Rotor-gene® Q	●	●		●			●
Thermo Scientific: PikoReal™	●	●		●			●
Illumina: The Eco™	●	●		●			●
Roche Applied Science: LightCycler® 480, LightCycler® Nano, LightCycler® 96	●	●		●			●
Roche Applied Science: LightCycler® 1.x, 2.0					●		

FIREPol is a trademark of Solis BioDyne. Applied Biosystems and StepOne are trademarks of Applied Biosystems LLC. QuantStudio, ViiA and Mustang Purple are trademarks of Life Technologies Corporation. Stratagene, Mx3000P, Mx3005P and Mx4000 are trademarks of Agilent Technologies Inc. Mastercycler is a trademark of Eppendorf AG. Rotor-Gene is a trademark of Qiagen Group. Bio-Rad, CFX96, CFX384, iQ, MyiQ, Opticon 2, Chromo4, MiniOpticon are the trademarks of Bio-Rad Laboratories. LightCycler and TaqMan are trademarks of Roche Molecular Systems Inc. The Eco is a trademark of Illumina Inc. PikoReal is a trademark of Thermo Fisher Scientific Inc. EvaGreen is a trademark of Biotium Inc. 3D model of FIREPol DNA Polymerase by Aare Abroi (Tartu University) and Sebastien Langui (Solis BioDyne), picture on the cover by Remo Savisaar.

qPCR Mix Compatibility Table: Probe-based qPCR Mixes

qPCR Platforms	5x HOT FIREPol® Probe qPCR Mix Plus (ROX)	5x HOT FIREPol® Probe qPCR Mix Plus (no ROX)	5x HOT FIREPol® Probe qPCR Mix Plus (Capillary)	5x HOT FIREPol® Probe Universal qPCR Mix	5x HOT FIREPol® Multiplex qPCR Mix	5x HOT FIREPol® Multiplex qPCR Mix (ROX)	5x HOT FIREPol® Multiplex qPCR Mix (Purple)
Applied Biosystems: 5700, 7000, 7300, 7700, 7900HT, StepOne™, StepOnePlus™	●			●		●	
Applied Biosystems: 7500, ViiA™7, QuantStudio™ 3*, 5, 6 Flex, 7 Flex, 12K Flex	●			●		●	●
Agilent/Stratagene: Mx3000PT™, Mx3005PT™, Mx4000™	●			●		●	
BioRad: iQ™5, MyiQ™, Chromo4™, Opticon®2; MiniOpticon®, CFX96™, CFX384™		●		●	●		
Bio Molecular Systems: Mic		●		●	●		
Eppendorf: Mastercycler® ep Realplex		●		●	●		
Qiagen: Rotor-gene® 3000, Rotor-gene® 6000, Rotor-gene® Q		●		●	●		
Thermo Scientific: PikoReal™		●		●	●		
Illumina: The Eco™		●		●	●		
Roche Applied Science: LightCycler® 480, LightCycler® Nano, LightCycler® 96		●		●	●		
Roche Applied Science: LightCycler® 1.x, 2.0			●				

* 5x HOT FIREPol® Multiplex qPCR Mix (Purple) is not compatible with Applied Biosystems QuantStudio™ 3.

5x HOT FIREPol® SolisGreen™ qPCR Mix

Description

Our newest real-time qPCR master mix based on DNA binding dye technology. SolisGreen™ is a new and unique dye with higher fluorescence, which enables detection of lower target concentrations.

Benefits

- improved detection of **low target concentrations**
- **higher** fluorescence level
- reaction set-up and shipment **without dry ice**
- compatible on **most cyclers** (see table on page 5)

Did you know?

SolisGreen™ and EvaGreen® dyes are detected in the same channel as SYBR Green I. You don't have to change any detection settings in your qPCR cycler.

Researchers already trust SolisGreen™

Reference:

“Excellent product quality along with affordable prices and committed customer service: these are the reasons why Solis BioDyne is our strategic enzyme supplier since many years. Now we continue our partnership with the robust, stable and very sensitive SolisGreen mixes for our qPCR product platform.”


DAVIDE ROASCHIO

Scientist in Product Development

Loewe Biochemica GmbH, Germany


HIGHLY COMPETITIVE

Amplification of six tenfold dilutions of human GAPDH cDNA with 5x HOT FIREPol® SolisGreen™ qPCR Mix (red) and a qPCR mix from another vendor (green). Reactions were performed on an Applied BioSystems™ QuantStudio™ 6.



EXCELLENT SENSITIVITY

Amplification of human PPIA cDNA with 5x HOT FIREPol® SolisGreen™ qPCR Mix on an Applied BioSystems™ QuantStudio™ 6 shows excellent linearity across different DNA concentrations. The amount of cDNA per reaction ranges from 100 ng to 0.1 pg.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® SolisGreen™ qPCR Mix	08-46-0000S (free sample)	50	0.2
	08-46-00001	250	1
	08-46-00008	2000	8
	08-46-00020	5000	20

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

5x HOT FIREPol® EvaGreen® qPCR Supermix

Description

Precisely-optimized real time qPCR master mix based on DNA binding dye technology. This master mix has been developed to give highly specific and sensitive results.

Benefits

- **high sensitivity** with low DNA concentrations
- **reduced primer dimer** formation
- reaction set-up and shipment **without dry ice**
- one qPCR mix for **all cyclers** (except capillary)
- contains **dUTP** to prevent cross-contamination with UNG treatment

Researchers already trust Supermix

Reference:

“ I use the 5x HOT FIREPol EvaGreen qPCR supermix to determine the expression level of low expressed genes after single cell RT in drop. This qPCR is for me a validation step before sending for sequencing. There are a lot of PCR inhibitors in my sample after RT in drop and my target cDNAs are in a very little amount. Indeed, I am studying low expressed genes and I am collecting just few thousand cells. However, I got good results with this qPCR mix, better than what I can obtain with a high fidelity polymerase. The melt curve are clean and I can see the variation in expression in between activated or not activated cells, even though those differences are very fine and even at high Ct values (around 35 cycles). Another advantage of this mix is the fact that it is 5x concentrated. I am doing 384 wells plate qPCR and a 5x mix is for sure more interesting and is cost-saving in comparison with a 2x mix like other suppliers provide. My colleagues also use this qPCR mix. They were so happy of the results they obtained with this mix that they ask for the enzyme alone for their PCR. We are now all using this mix for our qPCR. ”



SOPHIE FOULON

PhD Student

ESCP, France

TRUSTWORTHY PERFORMANCE

Amplification plots of tenfold dilution series for the human GAPDH gene performed on Applied Biosystems™ ViiA™7 (upper graph) and Roche LightCycler® 480 (lower graph). The amount of DNA per reaction ranges from 0.01 to 1.0 ng. The results show high linear range and high efficiency across a wide range of DNA concentrations on different qPCR platforms.



Selected publications:

- Eskla, K.L., et al. Free Radical Biology and Medicine (2018).
- Hubackova, S. et al. Cell Death and Differentiation (2018).
- Rockenbach, M.F., et al. Tropical Plant Pathology (2018).
- Kohler, V., et al. Nucleic Acids Research (2018).
- Bednarek, K., et al. Am J Cancer Res (2018).
- Maksimov, V., et al. PLOS ONE (2018).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® EvaGreen® qPCR Supermix	08-36-0000S (free sample)	50	0.2
	08-36-00001	250	1
	08-36-00008	2000	8
	08-36-00020	5000	20

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

5x HOT FIREPol® EvaGreen® qPCR Mix Plus

Description

Cost-effective real time qPCR master mix based on DNA binding dye technology. Provides highly reliable and reproducible results.

Benefits

- high sensitivity and specificity
- excellent efficiency
- reaction set-up and shipment **without dry ice**
- master mixes for ROX, no ROX and capillary cyclers

Researchers already trust EvaGreen® Plus

Reference:

“ We have been using HOT FIREPol® EvaGreen® qPCR Mix Plus (ROX) for almost four years. It works perfectly well as the competition (SYBR Green®) with the advantage of been thermostable. Moreover, the quality of the results coupled with its low price, make it the right choice for our objectives. Overall, Solis BioDyne is a great company with high quality products with immediate delivery. ”

LUIS FELIX



University of Trás-os-Montes and Alto Douro, Portugal

Selected publications:

- Weltner, J., et al. Nature Communications (2018).
- Lagos, C., et al. Clinical Immunology (2018).
- Reinhart, R., et al. Cell Death and Differentiation (2018).
- Hammoudi, V., et al. PLOS Genetics (2018).
- Choi, D., et al. Scientific Reports (2018).
- Sarmiento-Villamil, J.L., et al. Molecular Plant Pathology (2018).

EXCELLENT SENSITIVITY AND SPECIFICITY

The amplification of a 98 bp fragment of GAPDH gene exhibits sensitive and efficient reaction curves (upper graph) with highly specific peak in melt curve analysis (lower graph) using HOT FIREPol® EvaGreen® qPCR Mix Plus (no ROX). Amplification was performed on human genomic DNA using Rotor-Gene® 6000 qPCR cycler following cycling protocols recommended by the supplier.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (ROX)	08-24-0000S (free sample)	50	0.2
	08-24-00001	250	1
	08-24-00008	2000	8
	08-24-00020	5000	20
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (no ROX)	08-25-0000S (free sample)	50	0.2
	08-25-00001	250	1
	08-25-00008	2000	8
	08-25-00020	5000	20
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (Capillary)	08-26-0000S (free sample)	50	0.2
	08-26-00001	250	1
	08-26-00008	2000	8
	08-26-00020	5000	20

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

5x HOT FIREPol® EvaGreen® HRM Mix

Description

Reliable and sensitive real time qPCR master mix for high resolution melt (HRM) analysis.

Benefits

- reaction set-up and shipment **without dry ice**
- mixes available for ROX and no ROX cyclers
- **sensitive** EvaGreen® dye allows detection of DNA sequence variations

Did you know?

High resolution melt analysis can be used for SNP genotyping, discovering mutations, screening for heterozygosity, analysing DNA methylation.

Researchers already trust EvaGreen® HRM Mix

Reference:



“ We appreciate your products because most of our research takes place in Kenya and it can be difficult to ship items that are temperature sensitive. The Hot FIREPol® EvaGreen® HRM Mix is very effective and the prices are reasonable. ”

DR JEREMY HERREN

International Centre of Insect Physiology and Ecology, Kenya

SENSITIVE HRM GENOTYPING

High Resolution Melt Analysis was used to genotype a C insertion in BRCA1 gene, a breast cancer susceptibility gene, with HOT FIREPol® EvaGreen® HRM Mix (two graphs below). Reactions were performed on Corbett Rotor-Gene® 6000. Green lines represent wildtypes without an insertion, red lines represent a C insertion and blue line represents a patient with unknown phenotype.



No.	Color	Name	Genotype	Confi. %
37	Blue	Unknown phenotype	5382 wt	99,18
40	Green	Wildtype 1	5382 wt	97,33
41	Green	Wildtype 2	5382 wt	100,00
42	Red	Mutation 1	5382 mut insC	100,00
43	Red	Mutation 2	5382 mut insC	97,47

Selected publications:

- Morinha, F., et al. Conservation Genetics Resources (2018).
- Negrisolo, S., et al. European Journal of Human Genetics (2018).
- Yu, X., et al. Theoretical and Applied Genetics (2018).
- Grzegorzewska, A.E, et al. Vaccine (2018).
- Behboudi, H., et al. International Immunopharmacology (2018).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® EvaGreen® HRM Mix (ROX)	08-33-0000S (free sample)	50	0.2
	08-33-00001	250	1
	08-33-00008	2000	8
	08-33-00020	5000	20
5x HOT FIREPol® EvaGreen® HRM Mix (no ROX)	08-31-0000S (free sample)	50	0.2
	08-31-00001	250	1
	08-31-00008	2000	8
	08-31-00020	5000	20

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

NEW

5x HOT FIREPol® Multiplex qPCR Mix

Description

Probe-based qPCR master mix that has been optimized for highly sensitive and accurate quantification of up to 4 targets in a single reaction. This master mix was developed for TaqMan® probes but is also suitable for other hydrolysis probe types.

Benefits


- analyze **1-4 targets** in 1 reaction
- high specificity and sensitivity
- reaction set-up and shipment **without dry ice**
- robust amplification of **GC-rich** targets
- contains **dUTP** to prevent cross-contamination when used in combination with UNG
- wide instrument compatibility (see table on page 6)

Tip!

Reduce reagent cost and reaction set-up time by detecting multiple targets in a single reaction.


EXCELLENT FOR 4-PLEX ASSAYS

5x HOT FIREPol® Multiplex qPCR Mix (Purple) was used in 4-plex qPCR amplification with 4 tenfold serial dilutions of human gDNA (gDNA concentration in a reaction ranges from 10 ng/μl to 0.01 ng/μl). Reactions were performed with Applied BioSystems™ QuantStudio™ 6 cyclers using Purple dye for normalization.



SAME LEVEL OF SENSITIVITY WITH MULTIPLEXING

5x HOT FIREPol® Multiplex qPCR Mix (Purple) was used in 4-plex or 1-plex qPCR amplification with 5 tenfold serial dilutions of human cDNA (cDNA concentration in a reaction ranges from 10 ng/μl to 0.001 ng/μl). Reactions were performed with Applied BioSystems™ QuantStudio™ 6 cycler using Purple dye for normalization. The results show virtually identical Ct values for the multiplex and singleplex reactions across a wide template concentration range.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl	SIZE in ml
5x HOT FIREPol® Multiplex qPCR Mix	08-01-0000S (free sample)	50	0.2
	08-01-00001	250	1
	08-01-00008	2000	8
	08-01-00020	5000	20
5x HOT FIREPol® Multiplex qPCR Mix (ROX)*	08-02-0000S (free sample)	50	0.2
	08-02-00001	250	1
	08-02-00008	2000	8
	08-02-00020	5000	20
5x HOT FIREPol® Multiplex qPCR Mix (Purple)*	08-03-0000S (free sample)	50	0.2
	08-03-00001	250	1
	08-03-00008	2000	8
	08-03-00020	5000	20

* See the passive reference dye and probe reporter dye compatibility table on page 22.

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

5x HOT FIREPol® Probe Universal qPCR Mix

Description

Precisely-optimized real time qPCR master mix for probe based assays. This master mix has been developed for TaqMan® probes but is also suitable for other hydrolysis probe types.

Benefits

- suitable for singleplex and **duplex** assays
- high specificity and sensitivity
- reaction set-up and shipment **without dry ice**
- one qPCR mix for **all cyclers** (except capillary)
- superior results with templates with up to **75% GC** content
- contains **dUTP** to prevent cross-contamination when used with UNG

Tip!

Probe based qPCR is recommended over a dye based approach when specificity is especially important.

Researchers already trust Probe Universal

Reference:

“ We use Solis products in all our research groups. The most used product is HOT FIREpol Probe Universal qPCR mix. The product is very efficient and economic, offering the best cost benefit of the market. The fact their products are stable for 30 days at room temperature is another fantastic feature. ”

LAÍS MOREIRA GRANATO PHD


Post-Doc

Centro de Citricultura Sylvio Moreira, Brazil

Supplied by Sinapse Biotecnologia


QPCR PERFORMANCE IN A DUPLEX REACTION:

Two genes from human gDNA were amplified in duplex reaction using HOT FIREPol® Probe Universal qPCR Mix. Excellent results were obtained from four 10x dilutions (starting from 10 ng/μl) from both genes. BAIP3 (blue) with GC-content 70.3% and efficiency 100% and GAPDH (yellow) with GC-content 56.1% and efficiency 98.4%. Reactions were performed on Applied Biosystems ViiA™ 7 Real-Time PCR System.



HIGHLY COMPETITIVE QPCR MIX:

Four 10x dilutions of 197 bp long fragment of B4G4 gene with GC-content 75.6% were amplified from human gDNA using 5x HOT FIREPol® Probe Universal qPCR Mix (blue) and qPCR Mix from another vendor (green). Reactions were performed on Applied Biosystems ViiA™ 7 Real-Time PCR System following cycling protocol recommended by each supplier.



Selected publications:

- Dolci, M., et al. Environmental Pollution (2018).
- Lettlova, S., et al. Cellular Physiology and Biochemistry (2018).
- Schiro, G., et al. Journal of Fungi (2018).
- Okino, C.H., et al. Molecular Biology Reports (2018).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl	SIZE in ml
5x HOT FIREPol® Probe Universal qPCR Mix	08-17-0000S (free sample)	50	0.2
	08-17-00001	250	1
	08-17-00008	2000	8
	08-17-00020	5000	20

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

5x HOT FIREPol® Probe qPCR Mix Plus

Description

Cost-effective real time qPCR master mix for probe based qPCR assays. This master mix has been developed for TaqMan® probes but is also suitable for other hydrolysis probe types.

Benefits

- suitable for singleplex and **duplex** assays
- high specificity and sensitivity
- reaction set-up and shipment **without dry ice**
- master mixes for ROX, no ROX and capillary cyclers

Researchers already trust Probe Plus

Reference:

“ 5x HOT FIREPol® Probe qPCR Mix Plus (ROX) has been extensively used in our lab on a routine basis. Throughout the years, we have successfully analysed thousands of samples using this reliable, optimised, consumer- and budget-friendly solution provided by Solis BioDyne team.

”

MARINA GRIGOROVA


PhD, Research Fellow of Human Molecular Genetics
University of Tartu, Estonia

Selected publications:

- Weissova, K., et al. Sleep Medicine (2018).
- Gryndler, M., et al. Mycorrhiza (2018).
- Cilia, G., et al. European Journal of Protistology (2018).
- Podralska, M., et al. BMC Cancer (2018).


HIGHLY COMPETITIVE

Three tenfold dilutions of 72 bp fragment of albumin gene were amplified from human genomic DNA using HOT FIREPol® Probe qPCR Mix Plus (red) and a qPCR mix from Company A (green). Reactions were performed on Applied Biosystems 7900HT Real-Time PCR System following cycling protocols recommended by the supplier.



QPCR PERFORMANCE IN A DUPLEX REACTION

Amplification of FAM labelled target SNAI1 (orange) and VIC labelled reference gene HPRT (yellow) was performed in a single reaction using HOT FIREPol® Probe qPCR Mix Plus. This multiplex qPCR was carried out on three tenfold dilutions of human placental cDNA on Applied Biosystems 7900HT Real-Time PCR System.



Send your sample request to solis@sbd.ee

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® Probe qPCR Mix Plus (ROX)	08-14-0000S (free sample)	50	0.2
	08-14-00001	250	1
	08-14-00008	2000	8
	08-14-00020	5000	20
5x HOT FIREPol® Probe qPCR Mix Plus (no ROX)	08-15-0000S (free sample)	50	0.2
	08-15-00001	250	1
	08-15-00008	2000	8
	08-15-00020	5000	20
5x HOT FIREPol® Probe qPCR Mix Plus (Capillary)	08-16-0000S (free sample)	50	0.2
	08-16-00001	250	1
	08-16-00008	2000	8
	08-16-00020	5000	20

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** [solis.biodyne](https://www.skype.com/en/contacts/biodyne), via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

Product Selection Guide: Standard PCR Enzymes and Master Mixes

	HOT START	FIDELITY VS TAQ	READY TO LOAD	AMPLIFICATION RANGE	CLONING TYPE	ROBUST ON GC	PAGE(S)
HOT FIREPol® DNA Polymerase	●	1x		5 kb	TA	**	15
5x HOT FIREPol® GC Master Mix	●	1x		5 kb	TA	***	16
5x HOT FIREPol® Multiplex Mix	●	1x		3 kb	TA	**	17
5x HOT FIREPol® Multiplex Mix Ready to Load	●	1x	●	3 kb	TA	**	17
5x HOT FIREPol® Blend Master Mix	●	5x		5 kb	Blunt/TA	*	18-19
5x HOT FIREPol® Blend Master Mix Ready to Load	●	5x	●	5 kb	Blunt/TA	*	18-19
FIREPol® DNA Polymerase		1x		5 kb	TA	**	20
5x FIREPol® Master Mix		1x		5 kb	TA	*	21
5x FIREPol® Master Mix Ready to Load		1x	●	5 kb	TA	*	21

HOT FIREPol® DNA Polymerase

Description

Chemically modified hot-start version of the thermostable Taq DNA polymerase FIREPol®. This enzyme is activated only after heat treatment which prevents any unspecific polymerase activity at lower temperatures during reaction set-up. HOT FIREPol® DNA polymerase is supplied with 2 reaction buffers, 25 mM MgCl₂ and an additive for difficult templates.

Benefits

- increased specificity and sensitivity
- reduced primer dimer formation
- reaction set-up and shipment without dry ice

Researchers already trust HOT FIREPol®

Reference:


“ We had some problems with the implementation of a protocol, we tried for a long time with different enzymes without any positive results. We tested the HOT FIREPol® and it was the perfect troubleshooting, besides the great technical support received from Solis BioDyne. ”

DR MARIA JOSE SUAREZ

CIHATA, University of Costa Rica

HIGHLY COMPETITIVE

Four genes from human gDNA were amplified in multiplex reaction using HOT FIREPol® DNA Polymerase (lane 1-2) and two other hot start enzymes from company A (lane 3-4) and company B (lane 5-6). HOT FIREPol® DNA Polymerase performed well with all four genes in both 10x dilutions.



Selected publications:

- Hammer, Q., et al. Nature Immunology (2018).
- Burgstaller, J.P., et al. Nature Communications (2018).
- Gupta, P.R., et al. Human Molecular Genetics (2018).
- Giehr, P., et al. Nucleic Acids Research (2018).
- Šebest, L., et al. American Journal of Physical Anthropology (2018).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE in U
HOT FIREPol® DNA Polymerase (5 U/μl)	01-02-0000S (free sample)	100
	01-02-00500	500
	01-02-01000	1000

Products and samples

can be ordered via e-mail: info@solisbiodyne.com, via skype: solis.biodyne, via phone: +372 740 9960, or via our e-shop: solisbiodyne.com

5x HOT FIREPol® GC Master Mix

Description

PCR master mix that has been developed for working with difficult GC-rich templates and DNA secondary structures. The master mix contains hot-start Taq polymerase HOT FIREPol®, MgCl₂, dNTPs and a special buffer for high yield GC-rich amplification.

Benefits

- excellent amplification with templates up to **79% GC** content
- suitable for templates up to 5 kb
- reaction set-up and shipment **without dry ice**
- vials of 100% DMSO and 25 mM MgCl₂ enable **flexibility** in reaction optimization

Did you know?

GC-rich templates need a special PCR master mix because GC-rich DNA is more difficult to amplify. It forms stable secondary structures that are more resistant to denaturing and cause unspecific amplification.

Researchers already trust GC Master Mix

Reference:


“ In our lab, GC Master mix gave excellent results with low-abundance, difficult-to-amplify targets. Afterwards, these PCR products were cloned into expression vectors and sequenced - and vast majority of sequences were intact. So the GC Master mix has low mutation rate and is a good cloning tool as well. ”

DR ILLAR PATA
IVEX Lab, Estonia

Publication: Maksimov, M., et al. Neuropsychologia (2015).


AMPLICONS OF VARIOUS GC-CONTENT

12 GC-rich genes were amplified from human gDNA using HOT FIREPol® GC Master Mix. Final concentration of DNA template and DMSO was 1 ng/μl and 10% respectively. The Master Mix performed well on templates with up to 79% GC content.



AMPLICONS OF VARIOUS LENGTHS FROM GC-RICH TEMPLATE

GC-rich fragments of various length from human gDNA B4GN4 gene were amplified with HOT FIREPol® GC Master Mix. Final concentration of DNA template and DMSO was 1 ng/μl and 10% respectively. The Master Mix performed well with fragments of up to 3000 bp in length.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl	SIZE in ml
5x HOT FIREPol® GC Master Mix	04-33-00S15 (free sample)	25	0.1
	04-33-00115	250	1
	04-33-02015	5000	20

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

5x HOT FIREPol® MultiPlex Mix & MultiPlex Mix Ready to Load

Description

Precisely-optimized PCR master mix for multiplex PCR assays. This master mix contains thermostable hot-start Taq polymerase HOT FIREPol®, MgCl₂, dNTPs and buffer. You just need to add template, primers and water.

Benefits

- analyze **multiple targets** per reaction
- increased sensitivity and yield
- prevents **primer dimer** formation
- reaction set-up and shipment **without dry ice**

Researchers already trust MultiPlex Mix

Reference:

“ Convinced with the performance and quality of the product in multiple applications: robust enzyme activity and reproducible results in single and highly multiplexed PCRs. A “must-have” in the laboratory. ”


DR. SERGEY YAKUSHEV

Head of the laboratory

Microsynth, Switzerland

SENSITIVE AND SPECIFIC RESULTS

Different genes from human gDNA were amplified in multiplex reactions using 5x HOT FIREPol® MultiPlex Mix. Amplicons ranging from 122 bp to 1340 bp show similar yield and excellent separation with simultaneous amplification of even 18 targets.



Did you know?

Products specifically developed for multiplex assays contain sufficient amount of reaction components for accurate amplification of all targets.

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® MultiPlex Mix with 10 mM MgCl ₂	04-34-00S20 (free sample) 04-34-00120 04-34-02020	25 250 5000	0.1 1 20
5x HOT FIREPol® MultiPlex Mix Ready To Load with 10 mM MgCl ₂	04-36-00S20 (free sample) 04-36-00120 04-36-02020	25 250 5000	0.1 1 20

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

5x HOT FIREPol® Blend Master Mix & Blend Master Mix Ready to Load

Description

Precisely-optimized PCR master mix for more demanding PCR assays. In addition to the hot-start Taq polymerase HOT FIREPol® this master mix contains a proofreading enzyme which offers enhanced performance.

Benefits

- increased yield, sensitivity and specificity
- up to **5x higher fidelity**
- suitable for templates up to 5 kb
- reduced **primer dimer** formation
- reaction set-up and shipment **without dry ice**

Did you know?

MgCl₂ is needed for polymerase activity. Low MgCl₂ concentration can result in no PCR product, however too high MgCl₂ concentration can increase unspecific amplification. The optimal concentration varies within experiments but usually 1.5-2 mM is sufficient.

Researchers already trust Blend

Reference:

“ Solis BioDyne has proven to be a great company that has customer-oriented services in molecular work. 5x HOT FIREPol Blend Master Mix is convenient to use, store and produces high quality results. I know what I am saying because I have used this product in 2011, 2014 and now 2017, because it outperforms same products from other companies. ”


EMMY CHEPKOECH

PhD student

University of Eldoret, Kenya

AMPLICONS OF VARIOUS LENGTH FROM DIFFERENT TEMPLATES

Lines 1-4 present an excellent amplification of fragments of various length from λ DNA. Lines 5 and 6 show two different amplicons amplified from mouse genomic DNA. All these reactions were carried out using 5x HOT FIREPol® Blend Master Mix Ready to Load with 7.5 mM MgCl₂.



Lane	Template	Amplicon Length
1	λ DNA	499 bp
2	λ DNA	1003 bp
3	λ DNA	1998 bp
4	λ DNA	4991 bp
5	Mouse genomic DNA	808 bp
6	Mouse genomic DNA	3838 bp

Selected publications:

- Bahram, M., et al. Nature (2018).
- Otsing, E., et al. Soil Biology and Biochemistry (2018).
- Mullett, M.S., et al. Forest Pathology (2018).
- Strobel, A., et al. Environmental Toxicology and Chemistry (2018).
- Lange, M.K., et al. Veterinary Parasitology (2018).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® Blend Master Mix with 7.5 mM MgCl ₂	04-27-00S15 (free sample) 04-27-00115 04-27-02015	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix with 10 mM MgCl ₂	04-27-00S20 (free sample) 04-27-00120 04-27-02020	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix with 12.5 mM MgCl ₂	04-27-00S25 (free sample) 04-27-00125 04-27-02025	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix with 15 mM MgCl ₂	04-27-00S30 (free sample) 04-27-00130 04-27-02030	25 250 5000	0.1 1 20

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl	SIZE in ml
5x HOT FIREPol® Blend Master Mix Ready to Load with 7.5 mM MgCl ₂	04-25-00S15 (free sample) 04-25-00115 04-25-02015	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix Ready to Load with 10 mM MgCl ₂	04-25-00S20 (free sample) 04-25-00120 04-25-02020	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix Ready to Load with 12.5 mM MgCl ₂	04-25-00S25 (free sample) 04-25-00125 04-25-02025	25 250 5000	0.1 1 20
5x HOT FIREPol® Blend Master Mix Ready to Load with 15 mM MgCl ₂	04-25-00S30 (free sample) 04-25-00130 04-25-02030	25 250 5000	0.1 1 20

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** [solis.biodyne](https://www.skype.com/en/contacts/biodyne), via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

FIREPol® DNA Polymerase

Description

Genetically modified thermostable Taq DNA polymerase that provides robust and reproducible results. FIREPol® DNA polymerase is supplied with 2 reaction buffers, 25 mM MgCl₂ and an additive for difficult templates.

Benefits

- robust amplification for **routine** applications
- suitable for templates up to 5 kb
- reaction set-up and shipment **without dry ice**

Researchers already trust FIREPol®

Reference:

“ I found that for FIREPol® DNA Polymerase the quality was comparable to similar products even though the price was much cheaper for the Solis product. Therefore, Solis BioDyne are head and shoulders above their competitors when it comes to value for money which is especially important given the funding situation in these straitened times. ”

DR. GARY LOUGHRAN

Research Fellow


School of Biochemistry and Cell Biology
University College Cork, Ireland

Selected publications:

- Garcia-Martin, A.B., et al. Veterinary Microbiology (2018).
- Bhutada, G., et al. BBA - Molecular and Cell Biology of Lipids (2018).
- Babben, S., et al. BMC Genomics (2018).
- Buglione, M., et al. Mammalian Biology (2018).
- Zanatta, D.T., et al. Biological Journal of the Linnean Society (2018).


MOUSE GENOMIC DNA

1200 bp fragment of Beta-synuclein gene was amplified from mouse genomic DNA using FIREPol® DNA Polymerase with two different buffers: B (lane 1-3) and BD (lane 4-6). Template DNA was used at three tenfold dilutions starting from 1 ng/μl. FIREPol® DNA Polymerase was used at 0.04 U/μl.



PLANT GENOMIC DNA

672 bp fragment was amplified from barley genomic DNA using FIREPol® DNA Polymerase with two buffers: B (lane 1-3) and BD (lane 4-6). Template DNA was used at three tenfold dilutions starting from 1 ng/μl. The enzyme performed well even at a template concentration as low as 0.01 ng/μl. FIREPol® DNA Polymerase was used at 0.04 U/μl.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE in U
FIREPol® DNA Polymerase (5 U/μl)	01-01-0000S (free sample)	100
	01-01-00500	500
	01-01-01000	1000
	01-01-02000	2000

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

5x FIREPol® Master Mix & Master Mix Ready to Load

Description

Optimized ready-to-use PCR master mix for routine PCR assays. This master mix contains thermostable Taq polymerase FIREPol®, MgCl₂, dNTPs and buffer with detergent. You just need to add template, primers and water.

Benefits

- **ready to load** version allows you to load the PCR product on your gel straight after cycling
- suitable for templates up to 5 kb
- reaction set-up and shipment **without dry ice**
- all-in-one master mix format reduces pipetting errors and **saves time**

Researchers already trust FIREPol® Master Mix

Reference:

“ We are working with FIREPol® Master Mix for few years now and we are very satisfied. We are using it for genotyping (many different PCR) and we have excellent results. The ready to load format is very time saving. ”


CLAUDINE CORNELOUP

Engineer assistant

Platform for Experimental Biology on Mice,
Genotyping Laboratory, ENS de Lyon, France

PLANT GENOMIC DNA

672 bp fragment was amplified from barley genomic DNA using FIREPol® Master Mix (lane 1-3) and FIREPol® Master Mix Ready to Load (lane 4-6). Template DNA was used at three tenfold dilutions starting from 1 ng/μl. The Master Mixes performed well even at a template concentration as low as 0.01 ng/μl.



Selected publications:

- Tamm, M., et al. Ecological Indicators (2019).
- Fedele, S. et al. Scientific Reports (2017).
- Mahfouz, N., et al. Plos One (2017).
- Nummert, G., et al. Virology (2017).
- Sakmanoglu, A., et al. Journal of Microbiological Methods (2017).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl	SIZE in ml
5x FIREPol® Master Mix with 7.5 mM MgCl ₂	04-11-00S15 (free sample) 04-11-00115	25 250	0.1 1
5x FIREPol® Master Mix with 12.5 mM MgCl ₂	04-11-00S25 (free sample) 04-11-00125	25 250	0.1 1
5x FIREPol® Master Mix Ready to Load with 7.5 mM MgCl ₂	04-12-00S15 (free sample) 04-12-00115	25 250	0.1 1
5x FIREPol® Master Mix Ready to Load with 12.5 mM MgCl ₂	04-12-00S25 (free sample) 04-12-00125	25 250	0.1 1

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: [solis.biodyne](https://www.skype.com/en/contacts/biodyne), via **phone**: +372 740 9960, or via our **e-shop**: [solisbiodyne.com](https://www.solisbiodyne.com)

Product Selection Guide: One-step RT-PCR

		Nr of targets per reaction	GC-rich temp-lates	dUTP	Passive reference dye	Compatible cyclers	Incompatible probe reporter dyes	Page
Probe-based detection	SOLIScript™ 1-step Probe Kit	1-2	*	No	ROX	All cyclers	ROX JUN Texas Red	23
	SOLIScript™ 1-step Multiplex Probe Kit	1-4	***	Yes	None	All cyclers except Applied BioSystems™ and Agilent		24
	SOLIScript™ 1-step Multiplex Probe Kit (ROX)	1-4	***	Yes	ROX	Applied BioSystems™ and Agilent cyclers	ROX JUN Texas Red	24
	SOLIScript™ 1-step Multiplex Probe Kit (Purple)	1-4	***	Yes	Purple	Applied BioSystems™ cyclers with Mustang Purple® channel	Cy5	24
Dye-based detection	One-step RT-qPCR kit for dye-based detection*					–	–	–
End-point detection	One-step RT-PCR kit for end-point detection*					–	–	–

* Product description. These kits are currently in development and will be available soon.

NEW SOLIScript™ 1-step Probe Kit

Description

Convenient kit format for performing highly specific cDNA synthesis and probe based qPCR in a single tube. The kit contains qPCR mix, SOLIScript™ RT enzyme mix with RiboGrip™ RNase inhibitor, and nuclease free water.

Benefits

- cDNA synthesis up to 60°C for **superior specificity**
- suitable for singleplex and **duplex** assays
- reaction set-up and shipment **without dry ice**
- high specificity and sensitivity
- one kit for all cyclers

Tip!

Elevating reaction temperature enables highly specific primer annealing during reverse transcription.

Researchers already trust 1-step Probe Kit

Reference:

“ We needed to perform PCR in the point-of-use where cold storage is not available. SOLIScript 1-step Probe kit showed a good performance for detecting MS2 bacteriophage RNA after being stored at RT for 30 days. We have applied SOLIScript 1-step Probe kit to the detection of viruses in environmental samples obtaining satisfactory results. ”


DAVID AGUADO and SÍLVIA BOFILL-MAS

University of Barcelona, Spain

Supplied by Genycell Biotech

EXCELLENT QUANTIFICATION IN DUPLEX ASSAYS

SOLIScript™ 1-step Probe Kit was used to perform 2-plex one-step RT-qPCR with five tenfold serial dilutions of human total RNA (RNA amount ranges from 4000 pg/μl to 0.4 pg/μl per reaction). Reactions were performed with Applied BioSystems™ QuantStudio™ 6 cycler using ROX dye for normalization.



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl
SOLIScript™ 1-step Probe Kit	08-57-0000S (free sample) 08-57-00250	50 250

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

NEW SOLIScript™ 1-step Multiplex Probe Kit

Description


Convenient kit format for performing highly specific cDNA synthesis and probe based qPCR in a single tube. This kit was developed for sensitive and accurate RNA quantification with up to 4 targets in the same reaction. Multiplex compatibility and increased stability at room temperature make it ideal for high throughput RNA analysis. The kit contains qPCR mix optimized for multiplex reactions, SOLIScript™ RT enzyme mix with RiboGrip™ RNase inhibitor, and nuclease free water.

Benefits

- cDNA synthesis up to 60°C for **superior specificity**
- analyze **1-4 targets** in 1 reaction
- reaction set-up and shipment **without dry ice**
- robust amplification of **GC-rich** targets
- contains **dUTP** to prevent cross-contamination when used in combination with UNG
- wide instrument compatibility (see table on page 22)

EXCELLENT FOR 4-PLEX ASSAYS

SOLIScript™ 1-step Multiplex Probe Kit was used to perform 4-plex one-step RT-qPCR with five tenfold serial dilutions of human total RNA (RNA amount ranges from 4000 pg/μl to 0.4 pg/μl per reaction). Reactions were performed with Applied BioSystems™ QuantStudio™ 6 cyclor using Purple dye for normalization.



Did you know?

Products specifically developed for multiplex assays contain sufficient amount of reaction components for accurate amplification of all targets.

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 μl
SOLIScript™ 1-step Multiplex Probe Kit	08-55-0000S (free sample) 08-55-00250	50 250
SOLIScript™ 1-step Multiplex Probe Kit (ROX)	08-59-0000S (free sample) 08-59-00250	50 250
SOLIScript™ 1-step Multiplex Probe Kit (Purple)	08-61-0000S (free sample) 08-61-00250	50 250

Products and samples

can be ordered via **e-mail**: info@solisbiodyne.com, via **skype**: solis.biodyne, via **phone**: +372 740 9960, or via our **e-shop**: solisbiodyne.com

Product Selection Guide: First Strand cDNA synthesis

	SPECIFICATION	PRODUCT FORMAT	REACTION TEMPERATURE	RNASE H ACTIVITY	PAGE(S)
FIREScript™ RT cDNA Synthesis Mix	Excellent enzyme for standard applications	<ul style="list-style-type: none"> 3-vial-kit format 4 priming options^a 	37°C-60°C	Full	26-27
FIREScript™ RT cDNA Synthesis Kit		<ul style="list-style-type: none"> 7-vial-kit format 4 priming options^a 	37°C-60°C	Full	26-27
FIREScript™ Kit		<ul style="list-style-type: none"> enzyme and reaction buffer only 	37°C-60°C	Full	26-27
SOLIScript™ RT cDNA Synthesis Mix	Superior specificity with complex templates and specific primers	<ul style="list-style-type: none"> 3-vial-kit format gene-specific primers only^b 	37°C-60°C	Reduced	28
SOLIScript™ RT cDNA Synthesis Kit		<ul style="list-style-type: none"> 6-vial-kit format 2 priming options^c 	37°C-60°C	Reduced	28
SOLIScript™ Kit		<ul style="list-style-type: none"> enzyme and reaction buffer only 	37°C-60°C	Reduced	28

^a Primers provided with the product: oligo dT, random, oligo dT and random combined. Gene-specific primers to be supplied by the user.

^b Gene-specific primers to be supplied by the user.

^c Primer provided with the product: oligo dT. SOLIScript™ RT is not suitable with random primers. Gene-specific primers to be supplied by the user.

FIREScript™ Reverse Transcriptase

Description

Genetically modified MMLV-based reverse transcriptase with increased thermostability for better performance. FIREScript™ RT is extremely stable and will remain fully active for up to one month at room temperature. This RT contains a functional RNase H domain which can increase the sensitivity of RT-qPCR.

Benefits

- high specificity and yield
- wide reaction temperature **37°C - 60°C**
- fast **15 min** reaction time
- reaction set-up and shipment **without dry ice**
- available in convenient mix format and flexible kit format

Did you know?

A higher reaction temperature during reverse transcription denatures RNA secondary structures, which results in higher yields of full length cDNA.

Researchers already trust FIREScript™

Reference:

“ I used the FIREScript RT cDNA Synthesis Kit. The reverse transcription was done with 1 µg RNA to be transcript, 5 µM random primers, 500 µM dNTPs (mix). Synthesis was done following recommended quick protocol. My results were very good and I will replace my current product with FIREScript in the future. ”

VALERIE


Research technician

University of Basel, Switzerland

Supplied by LucernaChem AG


HIGHLY COMPETITIVE ENZYME

cDNA was synthesized with five tenfold human RNA dilutions using FIREScript™ (blue) and cDNA synthesis kits from competitor A (green) and competitor S3 (orange). Downstream qPCR reactions were performed with 5x HOT FIREPol® qPCR Supermix using β2M primers on an Applied Biosystems™ QuantStudio™ 6.



EXCEPTIONAL STABILITY

Two tenfold RNA dilutions were reverse transcribed to cDNA using FIREScript™ RT that had been stored at -20°C (blue upper graph) and FIREScript™ RT that had been stored at room temperature for 4 weeks (red lower graph). Downstream qPCR reactions were performed using 5x HOT FIREPol EvaGreen Supermix. The results are equal for both storage conditions.



Selected publications:

- Wojciechowicz, T., et al. Int J Food Sci Tech (2018).
- Popēna, I., et al. Cell Communication and Signaling (2018).
- Aufschnaiter, A. et al. Frontiers in Molecular Neuroscience (2018).

Send your sample request to info@solisbiodyne.com		
PRODUCT	CAT. NO.	RXN/20 µl
FIREScript™ RT cDNA synthesis MIX*	06-16-0000S (free sample)	20
FIREScript™ RT cDNA synthesis MIX with Oligo (dT) and Random primers	06-20-00100 06-20-00500	100 500
FIREScript™ RT cDNA synthesis MIX with Oligo (dT) primer	06-18-00100 06-18-00500	100 500
FIREScript™ RT cDNA synthesis MIX with Random primers	06-19-00100 06-19-00500	100 500
FIREScript™ RT cDNA synthesis MIX without primers	06-17-00100 06-17-00500	100 500
FIREScript™ RT cDNA synthesis KIT	06-15-0000S (free sample)	20
	06-15-00050	50
	06-15-00200	200
FIREScript™ KIT	06-13-0000S (free sample)	20
	06-13-00050	50
	06-13-00200	200

* The sample includes all 4 priming options. Gene-specific primers to be supplied by the user.

RiboGrip™ RNase Inhibitor

Description

Novel *in silico*-engineered, protein-based RNase inhibitor which inhibits RNase A. RiboGrip™ is extremely stable and will remain fully active for up to one month at room temperature.

Benefits

- reaction temperature up to 50°C
- reaction set-up and shipment **without dry ice**
- included in FIREScript™/SOLIScript™ RT cDNA Synthesis Kit and Mix

Send your sample request to info@solisbiodyne.com		
PRODUCT	CAT. NO.	RXN/20 µl
RiboGrip™ RNase Inhibitor	06-22-00400 (free sample)	20
	06-22-01000	50
	06-22-04000	200

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** [solis.biodyne](https://www.skype.com/en/contacts/solis-biodyne), via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

SOLIScript™ Reverse Transcriptase

Description

Novel *in silico*-engineered, thermostable reverse transcriptase with reduced RNase H activity. SOLIScript™ RT is extremely stable and will remain fully active for up to one month at room temperature.

Benefits

- **superior specificity** with gene-specific primers
- wide reaction temperature **37–60°C**
- reaction set-up and shipment **without dry ice**
- available in convenient mix format and flexible kit format

Did you know?

Always use experimental data to optimize the reaction temperature for your gene-specific primers in cDNA synthesis to find the best balance between specificity and yield.

Researchers already trust SOLIScript™

Reference:


“ SOLIScript is definitely the best RT enzyme in its price category. It is a consistent and processive enzyme. In our experiments, we work with heavily modified RNA species, and yet SOLIScript showed excellent results. ”

KIRILL JEFIMOV


University of Bergen, Norway

SUPERIOR SPECIFICITY


cDNA synthesis reactions were performed using a β 2M gene-specific primer at 55°C following protocols recommended by the suppliers. Downstream qPCR reactions were performed with 5x HOT FIREPol qPCR Supermix on an Applied Biosystems™ QuantStudio™ 6. Reactions in which cDNA was synthesized using SOLIScript™ RT (blue) have a highly specific qPCR melt curve compared to Competitor S3 (orange) and Competitor R (green).



■ Solis BioDyne SOLIScript™



■ Competitor S3



■ Competitor R

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	RXN/20 µl
SOLIScript™ RT cDNA synthesis MIX	06-37-0000S (free sample)	20
	06-37-00100	100
	06-37-00500	500
SOLIScript™ RT cDNA synthesis KIT	06-35-0000S (free sample)	50
	06-35-00050	50
	06-35-00200	200
SOLIScript™ KIT	06-33-0000S (free sample)	50
	06-33-00050	50
	06-33-00200	200

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

TERMIPol® DNA Polymerase

Description

Thermostable DNA polymerase that has an increased efficiency for incorporating unconventional nucleotides such as ddNTPs, acyNTPs or fluorescent nucleotides. TERMIPol® DNA Polymerase is supplied with a reaction buffer and 100 mM MgCl₂.

Benefits

- high efficiency for incorporating **unconventional nucleotides**
- assay success rate of 99% in MALDI-TOF
- robust and reliable
- reaction set-up and shipment **without dry ice**

Researchers already trust TERMIPol®

Reference:

“ Our group is using the TERMIPol® already for 10 years for primer extension reactions with subsequent HPLC separation. Compared to similar products on the market TERMIPol® incorporates ddNTPs with high efficiency and low error rates. We highly recommend using this enzyme for SNP genotyping or bisulfite-based single CpG screening, as low as 1.25 U are sufficient per reaction. Since no detergents are used in storage and reaction buffers, primer extension reactions can be loaded unpurified on HPLC systems which saves time and costs. We are using this enzyme frequently and experienced TERMIPol® as robust and reliable enzyme offering highly efficient and reproducible results. ”

DR. SASCHA TIERLING

Universität des Saarlandes, Germany

Did you know?

The ability to incorporate unconventional nucleotides makes TERMIPol® suitable for primer extension, MassARRAY and MALDI-TOF mass spectrometry.

Selected publications:

- Bormann, F., et al. International Journal of Cancer (2018).
- Royo, J.L., et al. Molecular and Cellular Probes (2015).
- Gorokhova, S.G., et al. J. Clin Exp Cardiol (2014).
- Thorkildsen, L.T., et al. Gastroenterology Research and Practice (2013).
- Ilina, E.N., et al. Front Microbiol (2013).
- Tierling, S., et al. International Journal of Cancer (2012).

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE in U
TERMIPol® DNA Polymerase (5 U/μl)	01-03-0000S (free sample)	500
	01-03-00500	500
	01-03-02000	2000
HOT TERMIPol® DNA Polymerase (5 U/μl)	01-06-0000S (free sample)	500
	01-06-00500	500
	01-06-02000	2000

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** [solis.biodyne](https://www.skype.com/en/contacts/biodyne), via **phone:** +372 740 9960, or via our **e-shop:** [solisbiodyne.com](https://www.solisbiodyne.com)

dNTP Mix and Set

Description

Solis BioDyne's dNTPs are chemically synthesized and have 99% purity determined by HPLC. You can use our dNTPs for a wide range of molecular biology applications.

dNTP Set

Separate vials of dATP, dTTP, dGTP and dCTP at 100 mM concentration.

dNTP Mix

One solution of dATP, dTTP, dGTP and dCTP at 20 mM concentration each.

dUTP

dUTP is available in a separate vial with a concentration of 100 mM.


“ In 2005 we started to use the Solis BioDyne dNTP Set in our lab. Comparing the performance of Solis BioDyne dNTPs with two other suppliers in a mutation detection assay, we found similar or even higher FRET signals in our analysed samples. Since then, we use the Solis BioDyne dNTP Set in our lab in a wide range of DNA and RNA amplification techniques like end point PCR, mutation detection in FRET assays, qPCR, high resolution melting analysis etc. ”

JUERGEN SIEVERTSEN

Bernhard Nocht Institute for Tropical Medicine (BNITM),
Germany


dATP

Formula: $C_{10}H_{12}N_5O_{12}P_3$ (Anion)
Formula Weight: 487.18 (Anion)




dCTP

Formula: $C_9H_{12}N_3O_{13}P_3$ (Anion)
Formula Weight: 463.15 (Anion)




dGTP

Formula: $C_{10}H_{12}N_5O_{13}P_3$ (Anion)
Formula Weight: 503.18 (Anion)



dTTP

Formula: $C_{10}H_{12}N_2O_{14}P_3$ (Anion)
Formula Weight: 478.16 (Anion)



Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE
dNTP Set	02-21-0001S (free sample) 02-21-00100 02-21-00400	4x1 µmol / 4x0.01 ml 4x25 µmol / 4x0.25 ml 4x100 µmol / 4x1 ml
dNTP Mix	02-31-0001S (free sample) 02-31-00020 02-31-00100	0.8 µmol / 0.01 ml 20 µmol / 0.25 ml 100 µmol / 1.25 ml
dUTP	02-41-0000S (free sample) 02-41-00025	2.5 µmol / 0.025 ml 25 µmol / 0.25 ml

Products and samples


can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

100 bp DNA Ladder

1 kb DNA Ladder

Description


Solis BioDyne DNA ladders are convenient ready-to-use molecular weight markers for DNA fragment size determination on gel electrophoresis. The ladders are supplied in a loading buffer and are stable at ambient temperature. The 1 kb DNA Ladder contains 13 discrete DNA fragments ranging from 300 bp to 10,000 bp. The 100 bp DNA Ladder contains 13 discrete DNA fragments ranging from 100 bp to 3,000 bp.



6x DNA Loading Dye Buffers

Description

6x DNA Loading Dye Buffers are used to prepare DNA markers and samples for loading on agarose or polyacrylamide gels. The optimized solutions contain different mixtures of three dyes: Bromophenol Blue, Xylene Cyanol FF and Orange G for visual tracking of DNA migration during electrophoresis. 6x DNA Loading Dye Buffers containing Orange G are recommended for the analysis of small DNA molecules and have no DNA masking during gel exposure to UV light. 6x DNA Loading Dye Buffer Blue and Dye Buffer Double Blue make pipetting visually easy with its dark blue color.



LOADING DYE BUFFERS

In 1% agarose gel 1x TBE, Xylene Cyanol FF migrates along with ~3500 bp fragments. Bromophenol Blue migrates along with ~300 bp fragments and Orange G migrates along with ~40 bp fragments.

Lane	DNA Loading Dye Buffer
1	Blue
2	Double Blue
3	Orange and Blue
4	Orange

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE
100 bp DNA Ladder Ready to Load	07-11-0000S (free sample) 07-11-00050	1.5 µg / 0.015 ml 50 µg / 0.5 ml
1 kb DNA Ladder Ready to Load	07-12-0000S (free sample) 07-12-00050	1.5 µg / 0.015 ml 50 µg / 0.5 ml
6x DNA Loading Dye Buffer Blue	07-01-0000S (free sample) 07-01-00001 07-01-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Double Blue	07-02-0000S (free sample) 07-02-00001 07-02-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Orange and Blue	07-03-0000S (free sample) 07-03-00001 07-03-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Orange	07-04-0000S (free sample) 07-04-00001 07-04-00010	0.1 ml 1 ml 10 ml

Products and samples

can be ordered via e-mail: info@solisbiodyne.com, via skype: [solis.biodyne](https://www.skype.com/en/contacts/solisbiodyne), via phone: +372 740 9960, or via our e-shop: solisbiodyne.com

10x GC-rich Enhancer



Description

10x GC-rich Enhancer is used as PCR additive for difficult GC-rich templates. The optimized solution modifies melting behavior of nucleic acids and often enhances amplification of suboptimal PCR systems with high degree of secondary structures and GC-rich regions.

10x GC-rich Enhancer should be used at a defined working concentration (1x, 2x or 3x solution) and only if non-specific amplification occurs.

Applications

- additive for PCR reaction



25 mM MgCl₂

Description

Magnesium Chloride (MgCl₂) is an important component of PCR reactions. Concentration of MgCl₂ should be optimized according to reaction conditions (primer, template, dNTP, polymerase concentration).

Applications

- optimization of PCR, qPCR and RT-PCR reactions
- all other molecular biology techniques where MgCl₂ is needed



PCR Grade Water

Description

PCR Grade Water is deionized and autoclaved water suitable for use in all experiments that require nuclease-free water. PCR Grade Water is prepared without chemical additives and it is DNase, RNase and nuclease-free.

Applications

- PCR, qPCR and RT-PCR
- all other molecular biology techniques where pure water is needed

Send your sample request to info@solisbiodyne.com

PRODUCT	CAT. NO.	SIZE in ml
10x GC-rich Enhancer	05-16-0000S (free sample)	0.1
	05-16-00010	1
	05-16-00050	5
	05-16-00200	20
25 mM MgCl ₂	05-11-00025	2.5
	05-11-00050	5
	05-11-00200	20
PCR Grade Water	water-025	25
	water-100	100
	water-500	500

Products and samples

can be ordered via **e-mail:** info@solisbiodyne.com, via **skype:** solis.biodyne, via **phone:** +372 740 9960, or via our **e-shop:** solisbiodyne.com

Ordering

All Solis BioDyne products are shipped at ambient temperature, without using dry ice.

Our products can withstand ambient temperature up to 1 month without any loss of activity. However, regular storage at -20°C is recommended.

» Free samples

Solis BioDyne provides free samples of the entire product range enabling our clients to thoroughly test our products.

How to Order

Orders can be placed:

- via E-Shop: solisbiodyne.com
- by emailing: info@solisbiodyne.com
- with your account manager or local distributor
- via fax: +372 740 2079

Required Information

Following information is required while placing an order:

- shipping and invoice address
- contact person's name and phone number
- VAT number (EU only)
- product name and corresponding catalogue number

Shipping Cost

Depending on the order amount a shipping cost may be added to the invoice. Please contact us for shipping cost quotation.

Customer Care

We are committed to providing our customers excellent service. All inquiries will be responded to within 1 business day at most. All technical questions will be given high priority and our full attention.

Please contact us through Skype: [support.sbd](https://www.skype.com/join/support.sbd) or via e-mail: info@solisbiodyne.com

Customized solutions

This product catalogue contains standard products, tube sizes and kits. Please contact us for customized solutions.

Shipping

Unless agreed otherwise, all shipments abroad will be arranged via express courier service. Orders are confirmed generally within 1 business day (Monday to Friday, 8AM to 5PM, UTC+2) after reception. In most cases orders are shipped within 1 or 2 business days.

Delivery documents and other charges

For non-EU shipments, please inform us of the documents required for shipments to your country. Solis BioDyne is not liable for import duties and taxes or delays caused by the brokerage procedure or other third parties.

Payment Options

Solis BioDyne accepts payments by:

- wire transfer, based on invoice
- PayPal, based on invoice or for orders placed through e-shop
- credit card (VISA or Master Card) for orders placed through our e-shop

Checks are not accepted as a payment method.


Please see our full ordering conditions on solisbiodyne.com

International Presence

Solis BioDyne currently has clients in more than 110 countries. We make our reagents accessible globally either by supplying them directly or in some areas by relying on local distributors who share our high standards in service and technical support. Please see the distributors section on next page or contact us directly to find the most convenient way to order in your region.



We are based in Estonia, a member of the European Union and €-zone



Distributors

AFRICA

Gene Life Sciences

Algeria
T: +213550537381
E: contact@genelifesciences.com
www.genelifesciences.com

Univers Biomedical Sarl

Burkina Faso
T: +226 50 48 36 16/17
E: bagre@unibio.bf
www.unibio.bf

DELTA Trading & Development

Egypt
T: +20 2 3573 6679
E: mostafa@deltatd.com
www.deltatd.com

Truspec Biotechnology

Kenya
T: +254 702995818
E: info@truspectd.com
www.truspectd.com

ASIA

Biogenuix Medsystems Pvt. Ltd.

India
T: +91 11 2561 2008
E: contact@biogenuix.com
www.biogenuix.com

Biolinkk

India
T: +91 1145035753
E: customerservice@biolinkk.com
www.biolinkk.com

PT. Medquest Jaya Global

Indonesia
T: +62 21 39842961
www.medquest.co.id

Next Gene Scientific

Malaysia
T: +603 5882 8650
E: sales@nextgene.my
www.nextgene.my

Med Lab Services

Pakistan
T: +92 51 4419341
E: mls@medlab.com.pk
www.medlab.com.pk

Genomicbase

Republic of Korea
T: +82 2 2215 4925
E: info@genomicbase.com
www.genomicbase.com

TLC Scientific Solutions Pte. Ltd.

Singapore
T: +656464 0261
E: orders@tlc-scientific.com
www.tlc-scientific.com

Krypto Genetics (PVT) Ltd.

Sri Lanka
T: +94 771 567 457
E: Orders@KryptoGenetics.com.lk
www.KryptoGenetics.com.lk

Omics Biotechnology

Taiwan
T: +886 2 8698 226
E: info@omicsbio.com
www.omicsbio.com

Life Science AP Co., Ltd

Thailand
T: +66 81 8296282
E: info@lifescienceap.com
www.lifescienceap.com

EUROPE

Medibena

Austria
T: +43 1 9906 497
E: info@medibena.at
www.medibena.at

Laborimpex SA

Belgium
T: +32 2 345 99 94
E: order@laborimpex.be
www.laborimpex.be

Biosistemi Bosnia I Hercegovina

Bosnia and Herzegovina
T: +387 33 264 316
www.biosistemigrupa.com

Biosistemi Hrvatska

Croatia
T: +385 (0)1 3460 839
E: biosistemi@biosistemi.hr
www.biosistemigrupa.com

Baria

Czech Republic
T: +420 774 227 421
F: +420 244 911 2 28
E: info@baria.cz
www.baria.cz

TAG Copenhagen A/S

Denmark
T: +45 321 322 00
E: oligo@tagc.com
www.tagc.dk

Labnet

Finland
T: +358 20 741 31 70
F: +358 20 741 31 89
E: labnet@labnet.fi
www.labnet.fi

Dominique Dutscher SAS

France
T: 03 88 59 33 90
E: hotline@dutscher.com
www.dutscher.com

Bioinnotech

Greece
T: +30 210 51 41 746
www.bioinnotech.gr

Bio-Kasztel

Hungary
T: +36 1 381 0694
E: info@kasztel.hu
www.kasztel.hu

CARLO ERBA Reagents S.r.l.

Italy
T: +39 02 93991003
E: laura.poletti@dgroup.it
www.carloerbareagents.com/it

SIA Adrona

Latvia
T: +371 67551894
E: info@adrona.lv
www.adrona.lv

Biosistemi MK d.o.o.e.l

Macedonia
T: +389 2 3216 622
www.biosistemigrupa.com

Bio-Connect B.V.

The Netherlands
T: +31 26 3264 450
E: info@bio-connect.nl
www.bio-connect.nl

Cytogen Polska Sp. z o.o.

Poland
T: +48 42 6300 598
E: cytogen@cytogen.com.pl
www.cytogen.com.pl

Bioportugal, Lda

Portugal
T: +351 226004800
E: bioportugal@bioportugal.pt
www.bioportugal.pt

Bio Zyme SRL

Romania
T: +40 264 523281
E: contact@biozyme.ro
www.biozyme.ro

Vivogen Srbija

Serbia
T: +381 11 6304 424
E: vivogen@vivogen.rs
www.biosistemigrupa.com

Ecoli s.r.o.

Slovakia
T: +421 0264 789 336
E: ecolli@ecoli.sk
www.ecoli.sk

Genycell Biotech

Spain
T: +34 902 194353
E: sales@genycell.com
www.genycell.com

LucernaChem AG

Switzerland
T: +41 41 4209 636
E: lucerna-chem@lucerna-chem.ch
www.lucerna-chem.ch

BM Yazılım Danis. ve

Laboratory Systems Ltd. Sti.
Turkey
T: +90 312 4472 280
E: bm@bmlabosis.com
www.bmlabosis.com

Newmarket Scientific

United Kingdom
T: +44 (0)1638 551500
E: tech@nktscientific.com
www.newmarketscientific.com

Tespro LLC

Ukraine
T: +380 44 220 110 5
E: tespro1@tespro.ua
www.tespro.com.ua

OCEANIA

Integrated Sciences Pty. Ltd.

Australia
T: +61 2 9417 7866
E: tech@integratedsci.com.au
www.integratedsci.com.au

dnature

New Zealand
T (toll-free): 0800 362 8873
E: info@dnature.co.nz
www.dnature.co.nz

MIDDLE EAST

BioConsult

Israel
T: +972 2 5667 043
E: leurer@bioconsult.co.il
www.bioconsult.co.il

Pishgam Biotech Company

Iran
T: +98 21 8801 4393
E: info@pishgambc.com
www.pishgambc.com

Al Genome Medical Company

Jordan
T: +962 652 33 670
E: info@cgenomix.com
www.cgenomix.com

Al Genome International Scientific

and Laboratory Products
Kingdom of Saudi Arabia
United Arab Emirates
T: +971557011398
E: info.uae@cgenomix.com

NORTH AMERICA

Mango Biotechnology LLC

USA
T: +1 650 5758 657
E: cs@mangobio.com
www.mangobio.com

SOUTH AMERICA

Biocientifica S.A.

Argentina
T: (54-11) 4857-5005
E: ventas@biocientifica.com.ar
www.biocientifica.com.ar

Sinapse Biotecnologia Ltda

Brazil
T: +55-11-2605-5655
E: sinapse@sinapsebiotecnologia.com.br
www.sinapsebiotecnologia.com.br

Fermelo Biotec

Chile
T/F: +562 2247 2978
E: info@fermelo.cl
www.fermelo.cl

SurGenoma S.A.S.

Colombia
T: +57 1 3441325
E: info@surgenoma.com
www.surgenoma.com



The most up-to-date list of distributors is available on our website solisbiodyne.com

Product List

qPCR Mixes		
	CAT. NO.	SIZE
5x HOT FIREPol® SolisGreen™ qPCR Mix	08-46-0000S 08-46-00001 08-46-00008 08-46-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® qPCR Supermix	08-36-0000S 08-36-00001 08-36-00008 08-36-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (ROX)	08-24-0000S 08-24-00001 08-24-00008 08-24-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (no ROX)	08-25-0000S 08-25-00001 08-25-00008 08-25-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® qPCR Mix Plus (Capillary)	08-26-0000S 08-26-00001 08-26-00008 08-26-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® HRM Mix (ROX)	08-33-0000S 08-33-00001 08-33-00008 08-33-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® EvaGreen® HRM Mix (no ROX)	08-31-0000S 08-31-00001 08-31-00008 08-31-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Multiplex qPCR Mix	08-01-0000S 08-01-00001 08-01-00008 08-01-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Multiplex qPCR Mix (ROX)	08-02-0000S 08-02-00001 08-02-00008 08-02-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Multiplex qPCR Mix (Purple)	08-03-0000S 08-03-00001 08-03-00008 08-03-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Probe Universal qPCR Mix	08-17-0000S 08-17-00001 08-17-00008 08-17-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Probe qPCR Mix Plus (ROX)	08-14-0000S 08-14-00001 08-14-00008 08-14-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Probe qPCR Mix Plus (no ROX)	08-15-0000S 08-15-00001 08-15-00008 08-15-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl
5x HOT FIREPol® Probe qPCR Mix Plus (Capillary)	08-16-0000S 08-16-00001 08-16-00008 08-16-00020	50 rxn/20 µl 250 rxn/20 µl 2000 rxn/20 µl 5000 rxn/20 µl

Regular PCR Enzyme and Mixes

	CAT. NO.	SIZE
FIREPol® DNA Polymerase (5 U/μl)	01-01-0000S	100 U
	01-01-00500	500 U
	01-01-01000	1000 U
	01-01-02000	2000 U
5x FIREPol® Master Mix with 7,5 mM MgCl ₂	04-11-00S15	25 rxn/20 μl
	04-11-00115	250 rxn/20 μl
5x FIREPol® Master Mix with 12.5 mM MgCl ₂	04-11-00S25	25 rxn/20 μl
	04-11-00125	250 rxn/20 μl
5x FIREPol® Master Mix Ready to Load with 7,5 mM MgCl ₂	04-12-00S15	25 rxn/20 μl
	04-12-00115	250 rxn/20 μl
5x FIREPol® Master Mix Ready to Load with 12,5 mM MgCl ₂	04-12-00S25	25 rxn/20 μl
	04-12-00125	250 rxn/20 μl

Hot-start PCR Enzyme and Mixes

	CAT. NO.	SIZE
HOT FIREPol® DNA Polymerase (5U/μl)	01-02-0000S	100 U
	01-02-00500	500 U
	01-02-01000	1000 U
5x HOT FIREPol® GC Master Mix	04-33-00S15	25 rxn/20 μl
	04-33-00115	250 rxn/20 μl
	04-33-02015	5000 rxn/20 μl
5x HOT FIREPol® MultiPlex Mix with 10 mM MgCl ₂	04-34-00S20	25 rxn/20 μl
	04-34-00120	250 rxn/20 μl
	04-34-02020	5000 rxn/20 μl
5x HOT FIREPol® MultiPlex Mix Ready To Load with 10 mM MgCl ₂	04-36-00S20	25 rxn/20 μl
	04-36-00120	250 rxn/20 μl
	04-36-02020	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix with 7.5 mM MgCl ₂	04-27-00S15	25 rxn/20 μl
	04-27-00115	250 rxn/20 μl
	04-27-02015	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix with 10 mM MgCl ₂	04-27-00S20	25 rxn/20 μl
	04-27-00120	250 rxn/20 μl
	04-27-02020	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix with 12.5 mM MgCl ₂	04-27-00S25	25 rxn/20 μl
	04-27-00125	250 rxn/20 μl
	04-27-02025	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix with 15 mM MgCl ₂	04-27-00S30	25 rxn/20 μl
	04-27-00130	250 rxn/20 μl
	04-27-02030	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix Ready to Load with 7.5 mM MgCl ₂	04-25-00S15	25 rxn/20 μl
	04-25-00115	250 rxn/20 μl
	04-25-02015	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix Ready to Load with 10 mM MgCl ₂	04-25-00S20	25 rxn/20 μl
	04-25-00120	250 rxn/20 μl
	04-25-02020	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix Ready to Load with 12.5 mM MgCl ₂	04-25-00S25	25 rxn/20 μl
	04-25-00125	250 rxn/20 μl
	04-25-02025	5000 rxn/20 μl
5x HOT FIREPol® Blend Master Mix Ready to Load with 15 mM MgCl ₂	04-25-00S30	25 rxn/20 μl
	04-25-00130	250 rxn/20 μl
	04-25-02030	5000 rxn/20 μl

First Strand cDNA Synthesis

	CAT. NO.	SIZE
FIREScript™ RT cDNA synthesis MIX sample	06-16-0000S	20 rxn/20 µl
FIREScript™ RT cDNA synthesis MIX with Oligo (dT) and Random primers	06-20-00100 06-20-00500	100 rxn/20 µl 500 rxn/20 µl
FIREScript™ RT cDNA synthesis MIX with Oligo (dT) primer	06-18-00100 06-18-00500	100 rxn/20 µl 500 rxn/20 µl
FIREScript™ RT cDNA synthesis MIX with Random primers	06-19-00100 06-19-00500	100 rxn/20 µl 500 rxn/20 µl
FIREScript™ RT cDNA synthesis MIX without primers	06-17-00100 06-17-00500	100 rxn/20 µl 500 rxn/20 µl
FIREScript™ RT cDNA synthesis KIT	06-15-0000S 06-15-00050 06-15-00200	20 rxn/20 µl 50 rxn/20 µl 200 rxn/20 µl
FIREScript™ KIT	06-13-0000S 06-13-00050 06-13-00200	20 rxn/20 µl 50 rxn/20 µl 200 rxn/20 µl
SOLIScript™ RT cDNA synthesis MIX	06-37-0000S 06-37-00100 06-37-00500	20 rxn/20 µl 100 rxn/20 µl 500 rxn/20 µl
SOLIScript™ RT cDNA synthesis KIT	06-35-0000S 06-35-00050 06-35-00200	50 rxn/20 µl 50 rxn/20 µl 200 rxn/20 µl
SOLIScript™ KIT	06-33-0000S 06-33-00050 06-33-00200	50 rxn/20 µl 50 rxn/20 µl 200 rxn/20 µl
RiboGrip™ RNase Inhibitor	06-22-00400 06-22-01000 06-22-04000	20 rxn/20 µl 50 rxn/20 µl 200 rxn/20 µl

One-step RT-qPCR Kits

	CAT. NO.	SIZE
SOLIScript™ 1-step Probe Kit	08-57-0000S 08-57-00250	50/20 µl 250/20 µl
SOLIScript™ 1-step Multiplex Probe Kit	08-55-0000S 08-55-00250	50/20 µl 250/20 µl
SOLIScript™ 1-step Multiplex Probe Kit (ROX)	08-59-0000S 08-59-00250	50/20 µl 250/20 µl
SOLIScript™ 1-step Multiplex Probe Kit (Purple)	08-61-0000S 08-61-00250	50/20 µl 250/20 µl

Other PCR Reagents		
	CAT. NO.	SIZE
TERMIPol® DNA Polymerase (5 U/μl)	01-03-0000S 01-03-00500 01-03-02000	500 U 500 U 2000 U
HOT TERMIPol® DNA Polymerase (5 U/μl)	01-06-0000S 01-06-00500 01-06-02000	500 U 500 U 2000 U
dNTP Set	02-21-0001S 02-21-00100 02-21-00400	4x1 μmol 4x25 μmol 4x100 μmol
dNTP Mix	02-31-0001S 02-31-00020 02-31-00100	0.8 μmol 20 μmol 100 μmol
dUTP	02-41-0000S 02-41-00025	2.5 μmol 25 μmol
100 bp DNA Ladder <i>Ready to Load</i>	07-11-0000S 07-11-00050	1.5 μg 50 μg
1 kb DNA Ladder <i>Ready to Load</i>	07-12-0000S 07-12-00050	1.5 μg 50 μg
6x DNA Loading Dye Buffer Blue	07-01-0000S 07-01-00001 07-01-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Double Blue	07-02-0000S 07-02-00001 07-02-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Orange and Blue	07-03-0000S 07-03-00001 07-03-00010	0.1 ml 1 ml 10 ml
6x DNA Loading Dye Buffer Orange	07-04-0000S 07-04-00001 07-04-00010	0.1 ml 1 ml 10 ml
10x GC-rich Enhancer	05-16-0000S 05-16-00010 05-16-00050 05-16-00200	0.1 ml 1 ml 5 ml 20 ml
25 mM MgCl ₂	05-11-00025 05-11-00050 05-11-00200	2.5 ml 5 ml 20 ml
PCR Grade Water	water-025 water-100 water-500	25 ml 100 ml 500 ml

SOLIS BIODYNE OÜ

A: Teaduspargi 9, 50411 Tartu, Estonia

REG NO: 10242922

VAT NO: EE100587614

E: info@solisbiodyne.com

T: +372 7409 960

F: +372 7402 079

SKYPE: solis.biodyne

B: Swedbank AS

IBAN code: EE692200221005142234

SWIFT: HABAEEX

Liivalaia 8, 15040 Tallinn, Estonia



European Union
European Regional
Development Fund



Investing
in your future